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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety precautions

EN

In order to prevent any possible risk of electric shock, fire and personal injury when the unit is being used, basic safety measures must always be taken, including:

- The heating plate can reach a temperature of up to 550°C. Care must be taken during the heating phase and also when the appliance is switched off as it takes time to cool down.
- The magnetic fields are active on magnetic and metallic parts (magnetic supports, CDs, floppy discs, cardiac stimulators, magnetic cards).
Keep these items away from the magnetic stirrer plate and the magnetic stir bars.
- The stirrer must not be used:
 - in hazard areas
 - in bain-marie
 - to stir combustible liquids that have a low combustion temperature.
- The vessel must be smaller in diameter than the ceramic plate to prevent any risk of tipping over.
- If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature – the use of Pyrex vessels is recommended.
- The following precautions must be taken before using the instrument:
 - position the instrument on a level surface.
 - position the instrument at least 30 cm from the wall
 - an operator must be present at all times when the instrument is in use
- The liquid must not come into contact with the electric power cable or with the electrical parts inside the instrument.
- Check that the power supply corresponds to the rating plate on the rear of the unit and that it is a good distance from the heating plate.
- Only correctly earthed three-terminal power cables must be used.
- Do not use the unit if it is not working correctly. In case of malfunctioning, contact your nearest assistance centre.
- Personal protective equipment must be compatible with the possible risks posed by the material being processed and the glass parts.
- Follow the cleaning instructions described in this manual

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1. Introduction

EN

The AREC heating magnetic stirrer with ceramic plate, is used in the laboratory for general use and for all those applications that require a precise regulation of the stirring speed and of the heating plate temperature.

The ceramic heating plate (1) means that higher temperatures can be obtained compared to the traditional aluminium heating plate and it takes less time to heat the sample. Ceramic is an inert and very hard material and is resistant to almost any type of chemical or mechanical aggression; its surface characteristics remain unaltered over time.

The AREC is used to mix liquids using a magnetically driven stir bar placed inside the vessel. The size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. It is therefore advisable to choose the most suitable magnetic stir bar in relation to the quantity and characteristics of the liquid and the type of vessel used. The stir bar which satisfies most applications is code A00001060 (10x40 mm).

The instrument is turned on using the ON/OFF button (2). The stirring speed (3) and the temperature (4) of the heating plate are adjusted using the relative knobs on the front panel, the values are indicated on the digital display (5).

NOTE: If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.


Fig. 1

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- AREC Heating Magnetic Stirrer
- Power cable
- Instruction manual

2.1 Electrical connection

After unpacking the instrument, place the unit on a laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply.

Connect the unit to the power supply using the cable supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

2.2 Start-up

After connecting the instrument to the power supply, place a flask containing the sample and a suitable magnetic stir bar on the heating plate.

2.3 Information about construction materials

Housing	Polymer
Front panel	PET
Feet	EPDM
Heating plate	Ceramic

3. Operating controls

Connect the unit to mains and turn it on using the on-off button. The green led indicates that the instrument is On. Start the stirring function by turning the stirrer knob "Stirrer rpm" on the front panel. Speeds of from 50 to 1500rpm can be selected using the analogical scale around the knob. Bright bars rotating on the display indicate that the instrument is running. The speed of rotation of these bars is proportional to the stirring speed.

NOTE: The rotating bars are not displayed when the temperature function is in use.

To start the heating function turn the heating knob "Heating °C" on the front panel. Temperatures of from 5 to 550°C can be selected at 5°C intervals. The temperature setting is shown on the display. When heating is turned Off (Heating knob OFF) and the heating plate temperature exceeds 50°C, the message "Hot" flashes on the display until the temperature of the heating plate falls to below 50°C (or until the instrument is turned off using the ON/OFF).

NOTE: This warning is not active if the instrument is not connected to the power supply.

The knobs on the front panel are easily accessible and are well away from the heating plate in order to ensure maximum operator safety as well as safeguarding the electronic components inside the unit.

Error messages

If an error message appears on the display, please contact VELP Scientifica's technical service department. When the display shows an error message, the stirring and heating functions stop automatically.

NOTE: To remove the error message disconnect the instrument from the power supply.

Message Cause
AL1 Thermocouple is out of range AL2 Excessive heating time AL3 The stirring system doesn't run correctly

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent.


5. Technical data

Power supply	F20500011 230V-50/60Hz ; F20510011 115V – 60Hz
Total electric power	800W
Dimensions (WxHxD)	203x93.5x344 mm (8x3.68x13.5 in)
Weight	3.5 Kg (7.7 lbs)
Speed range	50 - 1500 rpm
Speed Control	Electronic feed-back
Speed selection	Analogue
Overtemperature protection	Yes
Temperature range	5 - 550°C (41 - 1022°F)
Temperature control	Electronic probe
Temperature selection	5°C intervals
Max. load bearing capacity of plate	25kg
Ceramic plate dimensions	200x200mm
Construction material	Die-cast polymer
Environmental temp. range	5 - 40°C (41 - 104°F)
Storage temp. range	-10 - 60°C (14 - 140°F)
Max humidity	80%
Operating mode	Continuous
Electr. safety level CEI-EN60529	IP42
Pollution degree CEI-EN61010-1	2
Installation cat. CEI-EN61010-1	2
Max. altitude of installation	4000m (13123ft) above sea level
Fuse diam. 5x20mm	2x5A (2x8A – 115V/60Hz)

6. Accessories

A00001056	Magnetic stir bar Ø. 6x35mm
A00001060	Magnetic stir bar Ø 10x40mm*
A00001061	Magnetic stir bar Ø 9.5x60mm
A00001071	Hemispheric bowl for 250ml flasks
A00001072	Hemispheric bowl for 500ml flasks
A00001073	Hemispheric bowl for 1000ml flasks
A00001069	Support rod

* Magnetic stir bar recommended


7. Spare Parts

10002097 Knob

10000239 Foot

8. Wiring Diagram

- | | |
|--------------------------------------|----------------------|
| 1. ON/OFF button | 2. Temperature probe |
| 3. Display and potentiometer board | 4. Counter board |
| 5. Electrical resistance | 6. Electric motor |
| 7. Fuses 5x20 3.15A (5A - 115V/60Hz) | |

9. Warranty

The unit is guaranteed against production defects for **25 months** from our invoice date. In accordance with this guarantee VELP SCIENTIFICA undertakes to repair any instruments resulting as faulty due to the quality of the materials used or poor workmanship. Faults arising due to inexpert handling/use or carelessness will not be replaced or repaired under warranty.

For more details please contact your Distributor.

Exclusions:

The warranty will be considered null and void for faults resulting from:

- inexperience and carelessness of the customer.
- repairs, maintenance or replacements carried out by unauthorized third parties.
- use of non-original spare parts.

10. Declaration of conformity

We the manufacturer VELP SCIENTIFICA s.r.l.
address Via Stazione, 16
20865 USMATE (MB)
Italy

under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1 (2001)
EN 61326-1 (2006)
2002/95/EC (RoHS)
2002/96/EC (WEEE)

and satisfies the essential requirements of the following directives:

Machinery directive 2006/42/EC

Low voltage directive 2006/95/EC

Electromagnetic compatibility directive 2004/108/EC

plus modifications and that the documents listed in annex I are available at VELP's offices as foreseen by the machinery directive.

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